

Automatic separator by 3 part / Syringe needle, Hub, Syylinder
Separator the syringe needle from Syringe
KHS-01



KHS-01 will be the last solution against the Self-injury by Syringe in Hospital and Protection against the infectious diseases.

KHS-01 is the indispensable item in the hospital, sanitarium



Features and strong points

- Fast application
- User can remove Syylinder. Syring needle and Hub will fall into the separated bottles.
- Ultra Violet lamp will sterilizer to Separator and Syringe needle and Hub.
- Disused Syringe needle and Hub will be stored antiseptic and safety bottles.
- User will change only the bottles for the maintainance.
- Small and light weight
- Easy moving and installation



Structure of Syringer

1. Syylinder
2. Hub
3. Syringer needle
4. Safety cap



KHS-01 is the perfect processing the medical waste and safety, fast removal

KHS-01 is easy and safety

Procedure and operating



1. Insert the syringer needle into the interposing fixture



2. Working by light sensor



3. Separate the needle, hub and sylinder by the machinery in Separator



4. Store needle and hub in glass bottles for recycling



5. Separate sylinder



6. Waste management company will take care of the filled isolated container



Instruction of the bottles which filled the disused needles

1. Display the open 1st bottle into the KHS-01.
2. Display the open 2nd bottle bisde the 1st bottle.
3. Close the lip when the bottle is full.

Suggested the small disinfection oil into the bottles for



Head quarter Busan University Miryang Campus Business Enterprise Center #408, Chunghak-Ri,
Samrangjin-Myun, Miryang-City, Kyungnam, Korea, 627-706
Tel : 82-070-889-2448 / Fax : 82-55-356-2448

Factory 344-1, Shinpyung-dong, Saha-gu, Busan, Korea
Tel : 82-51-206-6165 / Fax : 82-51-291-7335